

## Four Drawer Medical Equipment Carts

Helping hospitals transition from volume to value base care.

**MedViron®** manufactures hospital furniture and equipment that assists in infection prevention and may reduce the number of HAIs, may nurture patient satisfaction and which could influence HCAHPS scores, and provides a lower cost of ownership with reduced maintenance through enhanced durability.

### Standard Features:

- Each drawer has a removable drawer liner for ease of cleaning.
- All surfaces are thermofilmed with durable KYDEX®. KYDEX® stands-up to strong hospital cleaning chemicals.
- Two-locking and two non-locking 5" poly casters.
- Full extension, steel drawer slides.
- Also available in 5 standard woodgrain colors.



Medviron LLC  
844.633.8476  
<mailto:customerservice@medviron.com>  
[medviron.com](http://medviron.com)

9380 Pentatech Drive  
PO Box 286  
Zeeland, MI 49464

Manufactured in the U.S.A.

### Dimensions

Surface: 32" wide x 19" deep

Height: 38"

Base: 34" wide x 21" deep



Each drawer has a removable drawer liner with a steel drawer frame and steel drawer slides.