


39" & 46" Overbed Tables



- A. Top surface is thermofilmed with KYDEX®.
 - 4 x thicker than standard thermofilm giving it superior impact strength and durability
 - Certified GREENGUARD Gold for children and schools
 - Withstands a wide range of cleaners and disinfectants
 - Withstands nervous disassembly by the patient.
- B. Top has eased edge for patient comfort. Spill containment area retains 40oz liquid.
- C. Activation lever conveniently located for patient and caregiver use and used only for lowering table height.
- D. Superior column, tested for 10,000 cycles with no effect on performance. Easy gas spring adjustment. Smooth column operation/height adjustment.

- E. MedViron Overbed Table with Shelves/Drawers aid hospitals looking to eliminate the bedside stand by providing additional storage that is easily accessible by the patient.

Having personal items conveniently within reach assists the patient's connection and control with their environment increasing patient satisfaction.
- F. The angle of the shelf is designed for easy retrieval of items so people with limited capabilities can slide objects on the shelf towards them and these objects naturally go into their hands.
- G. X-base designed for maximum stability has super low profile for 5th wheel beds. Dual casters provide smooth motion and effortless mobility.

Comprised of over 60% total recycled content and virtually 100% recyclable. Ships fully assembled in 100% recyclable packaging.

Infection Prevention:

All surfaces are easily accessible and easily cleaned. The table top is supported underneath with a plate which can be easily wiped down. MedViron Overbed Tables have seamless tops and storage areas. The shelves and drawers are removable for easy cleaning, making the interior of the storage pedestal completely accessible for environmental services staff.

39" & 46" Overbed Tables



Overbed Table - 39" Width,
Shown here in Summer Flame



Overbed Table - 2-shelf,
39" Width, Summer Flame,



Overbed Table - 1-shelf,
46" Width, Eiderdown



Overbed Table - drawer/shelf,
46" Width, Fusion Maple



Overbed Table - 2-drawer,
46" Width, Honey Maple



Overbed Table - 2-shelf,
46" Width, Fusion Maple,

Specifications:	46" width	39" width
Configurations:	single-shelf, 2-shelf, 2-drawer, drawer-shelf	no shelves, 2-shelves
Options:	cup holders, vanity mirror	cup holders, vanity mirror
Surface height:	29in to 46in (74cm to 1.2m)	29in to 46in (74cm to 1.2m)
Width:	46in (1.2m)	39in (99m)
Depth:	19in (48cm)	19in (48cm)
Weight:	84lbs (38kg)	57lbs (38kgs)
Drawer dimension:	9.5"W x 13.5"D x 3"H (24cm x 34cm x 8cm)	
Shelf dimension:	10.5"W x 16.5"D x 1.5"H (27cm x 42cm x 4cm)	4.5"W x 15.5"D x 1.5"H (12cm x 40cm x 4cm)
Warranty	Limited lifetime warranty.	Limited lifetime warranty.

Manufactured in the USA.