

32" Overbed Table



Infection Prevention: All surfaces of MedViron Overbed Tables are easily accessible and easily cleaned. MedViron Overbed Tables have seamless tops.

Top surface is thermofilmed with KYDEX® which is: 4x thicker than standard thermofilm giving it superior impact strength and durability, certified GREENGUARD® Gold for children and schools, withstands a wide range of cleaners and disinfectants, and withstands nervous disassembly by the patient.

The overbed table top has an eased edge for patient comfort. The spill containment area can retain 40oz of liquid.

The activation lever is conveniently located for patient and caregiver and is used only when lowering the table's height.

Dual casters provide smooth, quiet motion and nearly effortless mobility. The overbed table base has low a profile to fit under 5th wheel beds.

Can your hospital's OBT's do this? 315 pounds, evenly distributed.

MedViron's superior column offers smooth height adjustment. The column has been tested for 10,000 cycles with no effect on performance and is strong enough to support a 315 pound evenly distributed load.

MedViron overbed tables are comprised of over 60% recycled material and are virtually 100% recyclable.

Ships fully assembled in 100% recyclable packaging.

Manufactured in the USA.

