



HEIGHT ADJUSTABLE DELIVERY CART

The MedViron Height Adjustable Delivery Cart's ergonomic design facilitates infant care in hospital maternity environments.

INFECTION PREVENTION:

Seamless horizontal and vertical panels.

Removable proprietary KYDEX shelf makes cleaning simple.

DURABILITY:

Patented column offers superior strength, stability, and smooth height adjustment.

The column has been tested for 10,000 cycles with no effect on performance.

5-year gas spring warranty.

Horizontal surfaces are thermofilmed with KYDEX® giving them superior impact strength.

EASE OF USE:

Four height adjustment handles, one in each corner, can be used independently to raise or lower bassinet height.

Base panel design allows caregivers to get close to the infant.

2 locking and 2 non-locking casters provide stability and easy maneuverability.

Manufactured in the USA.



MDC-0000-0000



HEIGHT ADJUSTABLE DELIVERY CART



SURFACE MATERIAL: KYDEX® THERMOFILM

The Medviron delivery cart’s horizontal surfaces are thermofilmed with KYDEX® which is 4x thicker than standard thermofilm giving it superior impact strength and durability.

KYDEX is certified GREENGUARD Gold for Children and Schools®. KYDEX withstands a wide range of cleaners and disinfectants.

SPECIFICATIONS:

Width:	35.5in (90cm)
Drop leaf extended width:	49.5in (125cm)
Depth:	22in (56cm)
Surface low height	27in (18cm)
Surface high height	36in (91cm)

5 year gas spring warranty.
Limited lifetime warranty for other components.

STANDARD COLOR CHOICES:



Chocolate Pear




Fusion Maple



Honey Maple



Summer Flame



Eiderdown