



32" OVERBED TABLE



INFECTION PREVENTION:

All surfaces of Medviron Overbed Tables are easily accessible and easily cleaned. Medviron Overbed Tables have seamless tops with an eased edge for patient comfort.

DURABILITY:

Medviron's superior column offers smooth height adjustment. The column has been tested for 10,000 cycles with no effect on performance.

5-year gas spring warranty.

EASE OF USE:

The activation lever is conveniently located for patient and caregiver and is used only when lowering the table's height.

Dual casters provide smooth, quiet motion and nearly effortless mobility.

The base has low a profile to fit under 5th wheel beds.

The spill containment area can retain 40oz of liquid.

Ships fully assembled in 100% recyclable packaging.

OPTIONS:

cupholders, vanity mirror

Manufactured in the USA.

Can your facility's overbed tables do this?

315 pounds.

MedViron's patented column is strong enough to support 315lbs of evenly distributed weight without any effect on the performance of the table once the weight is removed.



STATEMENT OF LINE:



OBT32



OBT32 w/cupholders



OBT32 w/vanity mirror



OBT32 w/vanity & cupholders

SURFACE MATERIAL: KYDEX® THERMOFILM

The Medviron Overbed Table's top surface is thermofilmed with KYDEX® which is 4x thicker than standard thermofilm giving it superior impact strength and durability.

KYDEX is certified GREENGUARD Gold for Children and Schools®. KYDEX withstands a wide range of cleaners and disinfectants.

SPECIFICATIONS:

Width:	32" (80cm)
Depth:	19" (48cm)
Surface height:	27.5" to 45.5" (70cm to 1.1m)
Base height:	2.1" (5.3cm)
Weight:	51lbs (36kg)
Working weight capacity:	25lbs (11.3kg)
Maximum weight capacity:	100lbs (136kg)

**Limited lifetime warranty.
5-year warranty on gas spring.**

STANDARD COLOR CHOICES:



Chocolate Pear



Fusion Maple



Honey Maple



Summer Flame



Eiderdown