



32" OVERBED TABLE WITH SHELF



INFECTION PREVENTION:

MedViron Overbed Tables have seamless tops.

Shelf is removable for easy cleaning and makes the interior of the storage area completely accessible for environmental services.

The table top is supported underneath with a plate which can be easily wiped down.

DURABILITY:

MedViron's patented column offers superior strength, stability, and smooth height adjustment.

The column has been tested for 10,000 cycles with no effect on performance.

5-year gas spring warranty.

Top surface is thermofilmed with KYDEX® giving it superior impact strength.

Storage shelf made with KYDEX® giving it superior durability.

EASE OF USE:

The angle of the patented shelf is designed for easy retrieval of items.

The activation lever is conveniently located for patient and caregiver and is used only when lowering the table's height.

Dual casters provide smooth, quiet motion and nearly effortless mobility.

The base has low a profile to fit under 5th wheel beds.

The spill containment area can retain 40oz of liquid.

Ships fully assembled in 100% recyclable packaging.

Manufactured in the USA.

Can your facility's overbed tables do this?

315 pounds.

MedViron's patented column is strong enough to support 315lbs of evenly distributed weight without any effect on the performance of the table once the weight is removed.





32" OVERBED TABLE WITH SHELF



SURFACE MATERIAL: KYDEX® THERMOFILM

The Medviron Overbed Table's top surface is thermofilmed with KYDEX® which is 4x thicker than standard thermofilm giving it superior impact strength and durability.

KYDEX is certified GREENGUARD Gold for Children and Schools®. KYDEX withstands a wide range of cleaners and disinfectants.

SPECIFICATIONS:

Width:	32" (80cm)
Depth:	18" (42cm)
Surface height:	27.5" to 45.5" (70cm to 1.1m)
Base height:	2.1" (5.3cm)
Shelf dimensions:	24in W x 14in D x 1.125in H
Shelf opening	3in
Weight:	51lbs (36kg)
Working weight capacity:	25lbs (11.3kg)
Maximum weight capacity:	100lbs (136kg)
5 year gas spring warranty.	
Limited lifetime warranty for other components.	

STANDARD COLOR CHOICES:



Chocolate Pear




Fusion Maple



Honey Maple



Summer Flame



Eiderdown