



39" & 46" OVERBED TABLES

INFECTION PREVENTION:

All surfaces are easily accessible and easily cleaned. Seamless top and storage areas.

Shelves and drawers are removable for easy cleaning, making the interior of the storage pedestal completely accessible for environmental services.

DURABILITY:

Top surface is thermofilmed with KYDEX®:

KYDEX is 4x thicker than standard thermofilm giving it superior impact strength and durability. Certified GREENGUARD Gold for Children and Schools®. Withstands a wide range of cleaners and disinfectants.

MedViron's superior column has been tested for 10,000 cycles with no effect on performance.

EASE OF USE:

Top has eased edge for patient comfort. Spill containment area retains 40oz liquid.

Activation lever conveniently located for patient and caregiver use and used only for lowering table height.

MedViron Overbed Table with Shelves/Drawers aid hospitals looking to eliminate the bedside stand by providing additional storage that is easily accessible by the patient.

Having personal items conveniently within reach assists the patient's connection and control with their environment increasing patient satisfaction.

The angle of the shelf is designed for easy retrieval of items so people with limited capabilities can slide objects on the shelf towards them and these objects naturally go into their hands.

X-base designed for maximum stability has super low profile for 5th wheel beds. Dual casters provide smooth motion and effortless mobility.

Over 60% total recycled content and virtually 100% recyclable.

Ships fully assembled in 100% recyclable packaging.

Manufactured in the USA.



MEDVIRON 46" WIDE OVERBED TABLE



MEDVIRON 39" WIDE OVERBED TABLE

STATEMENT OF LINE:



39" Overbed Table



39" Overbed Table with 2 shelves



46" Overbed Table with 1-shelf



46" Overbed Table with 2 shelves



46" Overbed Table 1-drawer and 1-shelf



46" Overbed Table with 2-drawers

SPECIFICATIONS:	46in WIDTH	39in WIDTH
Configurations	Single-shelf, 2-shelves, 2-drawers, or drawer-shelf	No shelves, 2 shelves
Options	cup holders, vanity mirror	cup holders, vanity mirror
Surface Height	27.5in to 44.5in	27.5in to 44.5in
Width	46in (1.2m)	39in (.99m)
Depth	19in (.48m)	19in (.48m)
Ship Weight	84lbs (38kg)	57lbs (25kg)
Drawer Dimension	9.5in W x 13.5in D x 3in H (24cm x 34cm x 8cm)	
Shelf Dimension	10.5in W x 16.5in D x 1.5in H (27cm x 42m x 4cm)	4.5in W x 16.5in D x 1.5in H
Warranty	Lifetime limited warranty 5-year gas spring warranty	Lifetime limited warranty 5-year gas spring warranty