



## GUARDIAN DELIVERY CART



**Designed to accompany the Guardian Bassinet, the Guardian Delivery Cart offers the durability of KYDEX® thermoplastic on all horizontal surfaces.**

### **INFECTION PREVENTION:**

---

Seamless horizontal and vertical panels.

Large lower drawer has full extension drawer slides providing easy access for cleaning.

### **DURABILITY:**

---

The Guardian Delivery Cart offers the durability of KYDEX® thermoplastic on all surfaces, providing superior impact resistance.

### **EASE OF USE:**

---

The pass-through top drawer can be opened from front or rear of the cart.

Large lower drawer makes delivery packs easily accessible.

Dual drop-leaf top expands caregiver workspace when needed.

Almost effortless mobility on 5" casters – 2 locking and 2 non-locking.

The caster locks hold the cart firmly in place.

The curved base reduces the chances of caregivers hitting their shins while moving the cart.

**Manufactured in the USA.**



# GUARDIAN DELIVERY CART



### **SURFACE MATERIAL: KYDEX® THERMOFILM**

The Medvion delivery cart's horizontal surfaces are thermofilmed with KYDEX® which is 4x thicker than standard thermofilm giving it superior impact strength and durability.

KYDEX is certified GREENGUARD Gold for Children and Schools®. KYDEX withstands a wide range of cleaners and disinfectants.

### **SPECIFICATIONS:**

Surface Width:	29in (66cm)
Drop leaf extended width:	44in (1.1m)
Depth:	22in (61cm)
Surface height	36in (91cm)
Top drawer:	24in W x 18in L x 5in D
Bottom drawer:	24in W x 18in L x 13.5" D

Limited lifetime warranty.

### **STANDARD COLOR CHOICES:**



Chocolate Pear




Fusion Maple



Honey Maple



Summer Flame



Eiderdown