



The MedViron Fetal Monitor Work Station's design facilitates bedside documentation, fetal monitoring, and the storage and organization of equipment and supplies allowing clinicians to focus on the care of mothers and their babies.

The contemporary design and woodgrain finish provide a less "clinical" appearance, while giving the patient confidence that they are receiving competent care.

Dimensions can be customized. These carts can be made to fit the requirements of your staff and fetal monitoring system needs.

EASE OF USE:

Integrated pull-out keyboard tray for easy access.

2 locking and 2 non-locking 5" casters provide stability and easy maneuverability.

Ergonomic inset handles.

Ample storage space for supplies and accessories.

Drawers fully extend for easy access to stored items.

Manufactured in the USA.

SURFACE MATERIAL: KYDEX® THERMOFILM

The monitor cart's horizontal surfaces are thermofilmed with KYDEX® which is 4x thicker than standard thermofilm giving it superior impact strength and durability.

KYDEX is certified GREENGUARD Gold for Children and Schools®. KYDEX withstands a wide range of cleaners and disinfectants.

SPECIFICATIONS:

Surface	26in wide x 24in deep (66cm x 61cm)
Base	25in wide x 25in deep (63.5cm x 63.5cm)
Shelf	20.5in wide x 19.75in deep (52cm x 50cm)
Drawer Internal	17.75in wide x 16.5in long x 6in deep (45cm x 42cm x 15.24cm)
Pull out tray	19.25in wide x 18in deep (49cm x 46cm)
Surface height	48in high (122cm)
Casters	5in diameter (12.7cm)

Limited lifetime warranty.