

**INFECTION PREVENTION:**

MedViron Overbed Tables have seamless tops.

The table's lower shelf is removable for easy cleaning; the interior of the storage area is completely accessible for environmental services.

The table top is supported underneath with a plate which can be easily wiped down.

**DURABILITY:**

MedViron's superior column offers superior strength, stability and smooth height adjustment.

The column has been tested for 10,000 cycles with no effect on performance.

5-year gas spring warranty.

Top surface is thermofilmed with KYDEX® giving it superior impact strength.

Storage shelf is also made with KYDEX giving it superior durability.

**EASE OF USE:**

The angle of the patented shelf is designed for easy retrieval of items.

The activation lever is conveniently located for patient and caregiver and is used only when lowering the table's height.

Dual casters provide smooth, quiet motion and nearly effortless mobility.

The base has low a profile to fit under 5th wheel beds.

Top has an eased edge for patient comfort.

The spill containment area can retain 40oz of liquid.

Ships fully assembled in 100% recycled packaging.

**OPTIONS:**

Cupholders and vanity mirror

**Manufactured in the USA.**

**STATEMENT OF LINE:**



**OBT32 w/SHELF SPECIFICATIONS:**


Width:	32" (81 cm)
Depth:	18" (46 cm)
Surface height:	27.5" to 45.5" (70 cm to 1.1 m)
Base height:	2.1" (5.3 cm)
Weight:	51lbs (23 kg)
Shelf dimensions:	24" W x 14"D x 1.125"H
Shelf opening:	3" (7.5 cm)
Working weight capacity:	25lbs (11.3 kg)
Maximum weight capacity:	100lbs (45.4 kg)

Limited lifetime warranty.  
5-year warranty on gas spring.

**OBT32 SPECIFICATIONS:**

Width:	32" (81 cm)
Depth:	19" (48 cm)
Surface height:	27.5" to 45.5" (70 cm to 1.1m)
Base height:	2.1" (5.3 cm)
Weight:	45lbs (20.5 kg)
Working weight capacity:	25lbs (11.3 kg)
Maximum weight capacity:	100lbs (45.4 kg)

Limited lifetime warranty.  
5-year warranty on gas spring.



Chocolate Pear

Eiderdown

Fusion Maple

Honey Maple

Linen

Summer Flame