



FEATURES:

MedViron Medical Equipment Carts look like high quality furniture because they are high quality furniture. These modern carts are durable and offer an attractive alternative for your healthcare facility.

Available in 6 standard woodgrain or solid colors.

Swivel casters, 5in diameter, 2 locking and 2 non-locking.

Each drawer has a removable plastic drawer liner for ease of cleaning. These liners are interchangeable with the MedViron Guardian Bassinet.

Horizontal surfaces are thermofilmed with KYDEX®. KYDEX withstands degradation from strong hospital cleaning chemicals.

Ships fully assembled in 100% recyclable packaging.

Manufactured in the USA.

SPECIFICATIONS:

Top dimensions	34in x 19in
Height	38in
Base dimensions	34in x 21in
Drawer internal dimensions	26.5in x 13.25in x 6in
Hardware options	Black inset pulls or stainless steel bow pulls

Limited lifetime warranty.

Corporate Headquarters: 9380 Pentatech Drive, Zeeland, MI 49464

Website: medviron.com

Email: customerservice@medviron.com

Toll free: 844.633.8476