

## OVERBED TABLES



### Easy to lift, user-friendly overbed tables

- ✓ Single-touch vertical height adjustment, lift the table with just the touch of a finger
- ✓ Spill-containment area can retain up to 40oz of liquid
- ✓ Seamless top surface is thermofilmed with KYDEX® for superior impact strength and durability
- ✓ Conveniently located action lever makes it easy to lower the table height
- ✓ Dual casters for smooth motion
- ✓ Non-handed design, can be used on either side of the bed
- ✓ All MedViron Overbed Tables ship fully assembled in 100% recyclable packaging.
- ✓ Made in the USA.

## 32" OVERBED TABLE



## 39" OVERBED TABLE



## 46" OVERBED TABLE



### INFECTION PREVENTION:

All surfaces are easily accessible and easily cleaned. Seamless top and storage areas.

Shelves and drawers are removable for easy cleaning, making the interior of the storage pedestal completely accessible for environmental services.

### DURABILITY:

Top surface is thermofilmed with KYDEX®: KYDEX is 4x thicker than standard thermofilm giving it superior impact strength and durability.

KYDEX is certified GREENGUARD Gold for Children and Schools®. Withstands a wide range of cleaners and disinfectants.

MedViron's superior column has been tested for 10,000 cycles with no effect on performance.

### EASE OF USE:

Top has eased edge for patient comfort. Spill containment area retains 40oz liquid.

Activation lever conveniently located for patient and caregiver use and used only for lowering table height.

MedViron Overbed Table with Shelves/Drawers aid hospitals looking to eliminate the bedside stand by providing additional storage that is easily accessible by the patient.

Having personal items conveniently within reach assists the patient's connection and control with their environment increasing patient satisfaction.

The angle of the shelf is designed for easy retrieval of items so people with limited capabilities can slide objects on the shelf towards them and these objects naturally go into their hands.

X-base designed for maximum stability has super low profile for 5th wheel beds.

OPTIONAL FEATURES:



Dual cup holders



Vanity mirror



Drop-leaf top (46" top only)

PATIENT STORAGE OPTIONS:



32" OBT - Storage Shelf



39" OBT - 2 Shelves



46" OBT 1-drawer/1-shelf fully optioned with vanity and cup holders



46" OBT - Single Shelf



46" OBT - 2 Shelves



46" OBT - 2 Drawers



**SPECIFICATIONS:**

	<b>32" Overbed Tables</b>	<b>32" Overbed Tables w/shelf</b>
Options	cup holders, vanity mirror	cup holders, vanity mirror
Width	32" (80 cm)	32" (80 cm)
Depth	19" (46 cm)	18" (48 cm)
Surface height	27" to 43.25" (70cm to 1.1m)	27" to 43.25" (70 cm to 1.1m)
Base height	1.625" (5.3cm)	1.625" (5.3 cm)
Weight	45lbs (20.5 kg)	51lbs (23kg)
Shelf dimensions		24" W x 14"D x 1.125"H
Shelf opening		3" (7.5cm)
Working weight capacity	25lbs (11.3 kg)	25lbs (11.3 kg)
Maximum weight capacity	100lbs (45.4 kg)	100lbs (45.4 kg)

	<b>39" Overbed Tables</b>	<b>46" Overbed Tables</b>
Configurations	No shelves, 2 shelves	Single-shelf, 2-shelves, 2-drawers, or 1-drawer&1-shelf
Options	cup holders, vanity mirror	cup holders, vanity mirror
Width	39" (.99m)	46" (1.2m)
Depth	19in (.48m)	19" (.48m)
Surface Height	27.5" to 43.75"	27.5" to 43.75"
Base height	2.1" (5.3cm)	2.1" (5.3 cm)
Weight	57lbs (25kg)	84lbs (38kg)
Drawer Dimension		9.5in W x 13.5in D x 3in H (24cm x 34cm x 8cm)
Shelf Dimension	4.5" W x 16.5" D x 1.5" H	10.5" W x 16.5" D x 1.5" H (27cm x 42m x 4cm)
Warranty	Limited Lifetime warranty, 5-year gas spring warranty	

**Corporate Headquarters:**  
 9380 PentatechDrive, Zeeland, MI 49464  
**Website:** [medviron.com](http://medviron.com)  
**Email:** [customerservice@medviron.com](mailto:customerservice@medviron.com)  
**Toll free:** 844.633.8476

