



MEDVIRON HEIGHT-ADJUSTABLE MEDICAL CARTS

- ✓ Height-adjustable top surface
- ✓ Dual adjustment levers
- ✓ Modern design
- ✓ Woodgrain color options available
- ✓ Manufactured in the USA
- ✓ Casters by TENTE®



BETTER MOBILITY. BETTER LIFE.



FEATURES:

- The **MedViron Height-Adjustable Medical Cart** is immediately recognizable because of its unique height-adjustable upper work surface with integrated push handles combined with a stationary lower work surface with spill containment for two usable work areas on one cart. The top work surface can be raised from 32.75" to 48.625", while the lower work surface with spill containment is at 30.5".
- The cart has two work surface height-adjustment handles which can be used independently or simultaneously.
- Patented height-adjustment columns come with a 5 year gas spring warranty, an industry best.
- **The MedViron Height-Adjustable Medical Cart** offers the durability of KYDEX® thermoplastic on its horizontal surfaces, the drawer cabinet top and the cart base. The carts also have KYDEX thermoformed drawer fronts.
- The spill containment area of the secondary work surface retains 16oz of liquid.
- The carts base panel was designed to allow nursing staff to stand as close as possible to the work surface without hitting their shins on the cart.
- The **MedViron Height-Adjustable Medical Cart** has TENTE® casters: 2 swivel, 1 locking, and 1 steering/locking caster.
- Drawer box depths can be customized to customer's needs.
- Ships fully assembled in 100% recyclable packaging.
- Manufactured in the USA.

SPECIFICATIONS:

Surface Dimensions:

35.5" wide x 19.5" deep (90cm x 50cm)

Surface height - Low: 32.75" (76 cm)

Surface height - High: 48.625" (1.25 m)

Spill Containment Top Dimensions:

22" wide x 14.75" deep (56cm x 37cm)

Spill Containment Top Height:

30.5" (80 cm)

Base Dimensions:

36.2" wide x 21.7" wide (91cm x 55cm)



Corporate Headquarters:

9380 Pentatech Drive, Zeeland, MI 49464

Website: medviron.com

Email: customerservice@medviron.com

Toll free: 844.633.8476