



- ✓ Modern design
- ✓ Height adjustable
- ✓ Pass-through drawer
- ✓ Woodgrain color options available
- ✓ Manufactured in the USA
- ✓ Casters by TENTE®



ERGONOMIC DESIGN FEATURES:

- Four patented height-adjustment handles, one in each corner, to raise or lower the bassinet's height
- A pass-through drawer with removable drawer bin
- Stainless steel push handles
- A pass-thru / pull-out work surface that can be accessed from either side of the bassinet
- A durable PETG plastic tub with mattress that accommodates Trendelenburg and Reverse-Trendelenburg patient positioning



SURFACE MATERIAL: KYDEX® THERMOFILM

The MedViron **Rise II Bassinet** utilizes KYDEX® thermofilm, which is 4x thicker than standard thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD Gold for Children and Schools®. KYDEX withstands degradation from a wide range of harsh hospital cleaners and disinfectants.

DURABILITY:

MedViron's patented aluminum® columns offer superior strength, stability, and smooth height adjustment.

The pneumatic spring has a 5-year warranty.

The PETG plastic tub has 4 times greater impact resistance than acrylic tubs.

The bassinet utilizes KYDEX® thermofilm for exceptional impact strength. The pass-thru storage drawer is also made with KYDEX® giving it superior durability.

INFECTION PREVENTION:

A hospital's Environmental Services department is critical for infection prevention. MedViron's stainless steel push handle design makes the top surface of the bassinet completely accessible for environmental services staff to clean.

The bassinet's pass-thru drawer has a removable drawer tub with rounded inside corners for easy cleaning.

The **Rise II Bassinet** has seamless horizontal and vertical panels also aid in the ease of cleaning.

EASE OF USE:

Patented: Four height adjustment handles, one in each corner, can be used independently to raise or lower bassinet height.

The pass-thru pull-out work surface can be accessed from either side of the bassinet and has a groove to keep items from rolling off.

Base panel design allows caregivers to get close to the infant.

5" casters, 2 swivel, 1 locking, 1 steering and locking.

Ships fully assembled in recyclable packaging.

SPECIFICATIONS:

Overall Width	37.5"w
Top surface	36"w x 21.25" d
Surface height	30.5" low to 40.25" high
Base	35"w x 22"d
Drawer internal dimensions	22.5"w x 13.25"d x 6"h
Mattress dimensions	26" L x 12.5" W
Hardware options	Black inset pulls or stainless steel bow pulls
Limited lifetime warranty.	

