



## HIGHLIGHTS

- **Cleanable seamless surfaces**
- **Durable KYDEX**
- **Pull-out keyboard tray**
- **Solid and woodgrain color options**
- **Manufactured in the USA**

## ERGONOMIC DESIGN FEATURES

The Fetal Monitor Work Station’s design facilitates bedside documentation, fetal monitoring, and the storage and organization of equipment and supplies. The dimensions of the Fetal Monitor Work Station can be customized. The work station can be manufactured to fit the requirements of your staff and fetal monitoring system needs.





## SURFACE MATERIAL KYDEX® 3D LAMINATE

The Fetal Monitor Work Station's horizontal surfaces are thermoformed with KYDEX 3D laminate, which is 4x thicker than standard rigid thermofilm giving it superior impact strength and durability. KYDEX is certified GREENGUARD® Gold.

KYDEX withstands degradation from a wide range of hospital cleaners and disinfectants.

### DURABILITY

The Fetal Monitor Cart's horizontal surfaces are thermofilmed with KYDEX®, which is 4x thicker than standard thermofilm, giving it superior impact strength and durability. KYDEX withstands a wide range of cleaners and disinfectants.

### INFECTION PREVENTION EASE OF USE

The Fetal Monitor Cart has seamless horizontal and vertical panels that may aid in ease of cleaning.

The integrated pull-out keyboard tray is thermofilmed with KYDEX, making it durable and easily cleaned.

2 locking and 2 non-locking 5" casters provide stability and easy maneuverability.

Stainless steel bow drawer pulls.

Ample storage space for supplies and accessories. Drawers fully extend for easy access to stored items.

\*\*\*Monitors and monitor arms are not included.\*\*\*

Six standard colors available. Contact us for custom color needs.

### SPECIFICATIONS

Height	44.8in high
Base Dimensions	22.8in wide x 26.7in deep
Cabinet Width	22.2in
Casters	5in diameter
Limited Lifetime Warranty	Yes