



FEATURES:

The MedViron Medication & Supply Cart is immediately recognizable because of its unique Starboard top surface with a tempered glass insert on the left for viewing items inside the top left drawer. The top also features integrated push handles on both sides.

- Durable construction with an attractive selection of colors.
- Two top drawers with drawer fronts with locks, keyed alike. Choice of black inset pulls or stainless steel bow pulls.
- 4 smaller easy-to-clean removable drawers with 2 large removable drawer bins below.
- Swivel casters, 5in diameter, 2 locking and 2 non-locking.
- Base panel is thermoformed with KYDEX® for durability and impact resistance.
- Ships fully assembled in 100% recyclable packaging. Manufactured in the USA.

SPECIFICATIONS

Cart dimensions	32in wide x 21.5in deep x 40in high
Small drawer dimensions	10.5in wide x 15in long x 3in deep
Large drawer dimensions	23in wide x 15in long x 6in deep
Colors	Top: white Sides, back and base: standard wood grain or solid color

