



HIGHLIGHTS

- **Cleanable seamless surface**
- **Integrated push handles**
- **Spill containment groove**
- **Optional locking door**
- **No assembly required**
- **Manufactured in the USA**

DESIGN FEATURES

Combining durability and convenience, these rugged Courier Carts make it easy to transport specimens throughout your medical facility. Four 5" casters move smoothly and quietly over virtually any surface, and two routed handles make maneuverability a breeze.



MADE IN U.S.A.



SURFACE MATERIAL

Starboard®

The Courier Specimen Transport Cart is immediately recognizable because of its unique Starboard top work surface with integrated push handles on either side. This durable polymer surface offers high impact and chemical resistance, is easy to clean and contains no volatile organic compounds (VOCs).

DURABILITY

Our Courier Specimen Transport Cart's durable design stands up to use throughout your entire facility.

This durable polymer top surface offers high impact resistance and chemical resistance, is easy to clean, and contains no volatile organic compounds (VOCs).

INFECTION PREVENTION

A hospital's Environmental Services department is critical for infection prevention. The surfaces of the cart are completely accessible for environmental services staff to clean.

All the cart's drawer bins are removable for easy cleaning. Each drawer bin has rounded inside corners for ease of cleaning.

The cart has seamless horizontal and vertical panels that may also aid in the ease of cleaning.

EASE OF USE

Available with or without a locking door.

A perimeter groove on the work surface prevents potential spills and specimen roll-offs.

Drawers easily slide out for quick, convenient access to contents, and are excellent for storing specimen bags and totes, and tube racks.

Four 5" casters (two locking). Designed for smooth, quiet rolling.

Integrated handles on either side of the work surface make the cart simple to move.

Ships fully assembled in 100% recyclable packaging.

SPECIFICATIONS

| | |
|--------------------|-----------------------|
| Overall Dimensions | 25" Wide x 21.3" Deep |
|--------------------|-----------------------|

| | |
|----------------|-----|
| Surface Height | 37" |
|----------------|-----|