



HIGHLIGHTS

- **Cleanable seamless surfaces**
- **Durable KYDEX®**
- **Impact resistant PETG tub**
- **Removable drawer liners**
- **Solid and woodgrain color options**
- **Manufactured in the USA**

ERGONOMIC DESIGN FEATURES

Base panel reduces the chances of caregivers hitting their shins.

Full-extension drawer slides.

Push handle height enables proper pushing posture.

Patented tub design allows for easier access to patient.

Trendelenburg and reverse Trendelenburg positioning.





SURFACE MATERIAL: **KYDEX® 3D LAMINATE**

The Guardian Bassinet utilizes KYDEX 3D laminate on all horizontal surfaces, which is 4x thicker than standard thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD® Gold. KYDEX withstands degradation from a wide range of hospital cleaners and disinfectants.

DURABILITY

KYDEX, with its superior impact resistance and durability, is thermoformed to all horizontal surfaces, including top, base panel and chart board. The PETG plastic tub has 4 times greater impact resistance than acrylic. The corners are reinforced with additional material to reduce the likelihood of corner breakage. Guardian Bassinets have stainless steel push handles.

INFECTION PREVENTION

Seamless KYDEX surfaces withstand commonly used hospital cleaners and disinfectants. Stainless steel push handles eliminate the nooks and hard to get to areas to make all of the surface area accessible and free from obstructions, allowing easy cleaning. Removable drawer liners have rounded corners that are easily wiped down and cleaned.

EASE OF USE

Stainless steel push handles allow greater view of infant. Patented tub design allows for easier access to patient. Base panel designed to reduce the chances of caregivers hitting their shins. Groove in the chart board keeps item from rolling off. Trendelenburg positioning for medical procedures. Full-extension drawer slides.

SPECIFICATIONS

Surface Height	30in (76 cm)
Handle Height	36in (92 cm)
Cabinet Width	34in (86 cm)
Cabinet Depth	21.75in (53 cm)
Mattress Dimensions	12.5in Wide x 26in long
Weight	112lbs (51 kg)