



## HIGHLIGHTS

- **Cleanable seamless surfaces**
- **Durable KYDEX**
- **Trendelenburg and Reverse Trendelenburg positioning**
- **Pass-through drawers**
- **Solid and woodgrain color options**
- **Manufactured in the USA**

## ERGONOMIC DESIGN FEATURES

Height is adjustable to the appropriate stature of each user and recovering C-section mothers.

Patented columns facilitate smooth height adjustment.

Patented height adjustment handles, one in each corner can be used independently to raise or lower the bassinet's surface height.

Base panel design allows caregivers to stand close to the work surface and closer to the newborn.





## SURFACE MATERIAL KYDEX® 3D LAMINATE

MedViron’s Rise Bassinet utilizes KYDEX 3D laminate, which is 4x thicker than standard rigid thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD® Gold.

KYDEX withstands degradation from a wide range of hospital cleaners and disinfectants.

## DURABILITY

MedViron’s patented aluminum columns offer superior strength, stability, and smooth height adjustment.

The pneumatic spring has a 5-year warranty.

The PETG plastic tub has 4 times greater impact resistance than acrylic tubs.

The bassinet utilizes KYDEX® 3D laminate for exceptional impact strength. The pass-thru storage drawers are also made with KYDEX® giving them superior durability.

## INFECTION PREVENTION

A hospital’s Environmental Services department is critical for infection prevention. MedViron’s stainless steel push handle design makes the top surface of the bassinet completely accessible for environmental services staff to clean.

The bassinet’s pass-thru drawer has a removable drawer tub with rounded inside corners for easy cleaning. The Rise Bassinet has seamless horizontal and vertical panels that may aid in the ease of cleaning.

## EASE OF USE

Patented: Four height adjustment handles, one in each corner, can be used independently to raise or lower bassinet height.

The optional flip-up work surface has a spill containment area for retaining liquids.

Base panel design allows caregivers to get close to the infant.

5” casters: 1 swivel, 2 locking, 1 steering & locking.

Ships fully assembled in recyclable packaging.

## SPECIFICATIONS

|                            |                                       |
|----------------------------|---------------------------------------|
| Surface Dimensions         | 36in wide x 21.25 deep                |
| Surface Height             | 30.5in low to 40.25 high              |
| Overall Width              | 37.5in wide                           |
| Mattress Dimensions        | 26in long x 12.5in wide               |
| Drawer Internal Dimensions | 22.5in wide x 13.25in deep x 6in high |