



**HIGHLIGHTS**

- **Cleanable seamless surfaces**
- **Durable KYDEX**
- **Patented Height Adjustment**
- **Trendelenburg and Reverse Trendelenburg positioning**
- **Pass-through drawer and work surface**
- **Solid and woodgrain options available**
- **Manufactured in the USA**

**ERGONOMIC DESIGN FEATURES**

Height adjustable to the appropriate stature of each user and recovering C-section mothers.

Patented column facilitates smooth height adjustment.

Patented height adjustment handles, one in each corner, can be used independently to raise or lower the Bassinet’s surface height.

Base panel design allows caregivers to stand close to the work surface and closer to the newborn.





## SURFACE MATERIAL: KYDEX 3D LAMINATE

The medViron Rise II Bassinet utilizes KYDEX 3D laminate, which is 4x thicker than standard rigid thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD Gold. KYDEX withstands degradation from a wide range of hospital cleaners and disinfectants.

### DURABILITY

MedViron's patented aluminum columns offer superior strength, stability, and smooth height adjustment.

The pneumatic spring has a 5-year warranty.

The PETG plastic tub has 4 times greater impact resistance than acrylic tubs.

The bassinet utilizes KYDEX® 3D laminate for exceptional impact strength. The pass-thru storage drawer is also made with KYDEX® giving it superior durability.

### INFECTION PREVENTION

A hospital's Environmental Services department is critical for infection prevention. MedViron's stainless steel push handle design makes the top surface of the bassinet completely accessible for environmental services staff to clean.

The bassinet's pass-thru drawer has a removable drawer tub with rounded inside corners for easy cleaning.

The Rise II Bassinet has seamless horizontal and vertical panels also aid in the ease of cleaning.

### EASE OF USE

Patented: Four height adjustment handles, one in each corner, can be used independently to raise or lower bassinet height.

The pass-thru pull-out work surface can be accessed from either side of the bassinet and has a groove to keep items from rolling off.

Base panel design allows caregivers to get close to the infant.

5" casters, 2 swivel, 1 locking, 1 steering and locking.

Ships fully assembled in recyclable packaging.

### SPECIFICATIONS

<b>Overall Width</b>	37.5" W
<b>Top Surface</b>	36" W x 21.25" D
<b>Surface Height</b>	30.5" low to 40.25" high
<b>Base</b>	35" W x 22" D
<b>Drawer Internal Dimensions</b>	22.5" W x 13.25" D x 6" H
<b>Mattress Dimensions</b>	26" L x 12.5" W
<b>Hardware Options</b>	Black Inset Pulls or Stainless Steel Bow Pulls

**Limited Lifetime Warranty**

