



HIGHLIGHTS

- Cleanable seamless surfaces
- Durable KYDEX® top surface
- Patented height adjustment
- Limited lifetime warranty
- Manufactured in the USA

DESIGN FEATURES

Patented column facilitates smooth height adjustment.

Height adjustment allows user to adjust to fit the task.

Dual casters provide smooth motion and mobility.

Top has eased edge for user comfort.





SURFACE MATERIAL: KYDEX® 3D LAMINATE

The medViron® Roam & Surf work surfaces utilize KYDEX® 3D laminate on the top surfaces, which is 4x thicker than standard rigid thermofilm, giving it superior impact strength and durability. KYDEX® is certified GREENGUARD® Gold. KYDEX® withstands degradation from a wide range of hospital cleaners and disinfectants.

DURABILITY

MedViron's patented aluminum column offers superior strength, stability, and smooth height adjustment.

Height adjustment column has been tested for 10,000 cycles with no effect on performance

Limited lifetime warranty. Pneumatic spring has a 5-year warranty.



Roam work surfaces are available with an optional drawer

INFECTION PREVENTION

A hospital's Environmental Services department is critical for infection prevention. Surfaces that are easy to clean are more likely to be cleaned.

Seamless thermoformed KYDEX® withstands commonly used hospital cleaners and disinfectants.

Top surface is supported underneath with a plate which is free from obstructions and is easily wiped down.

Optional drawer is removable for easy cleaning .

EASE OF USE

Patented column facilitates smooth height adjustment.

Height adjustment allows user to adjust to fit the task.

Versatile, to be used by caregivers, patients, or family members.

Ships fully assembled in recyclable packaging.

Roam and Surf work surface tops are available in solid or woodgrain options.

Visit our website for all color options.



SPECIFICATIONS

	ROAM	SURF
Width	23.5" (60 cm)	35" (89 cm)
Depth	20" (51 cm)	19" (48 cm)
Surface Height	29" to 46" (74 cm to 1.2 m)	29" to 46" (74 cm to 1.2 m)
Options	Drawer, Drawer w/ lock	—
Warranty	Limited lifetime/5-year gas spring	Limited lifetime/5-year gas spring