



## HIGHLIGHTS

- **Cleanable seamless surfaces**
- **Durable KYDEX**
- **Dual drop-leaf top surface**
- **Pass-through drawer**
- **Solid and woodgrain color options**
- **Manufactured in the USA**

## ERGONOMIC DESIGN FEATURES

Pass-through drawer can be opened from the front or rear of the cart.

A large drawer makes delivery packs easily accessible.

Dual drop-leaf top expands caregiver workspace when needed.

Curved base reduces chances of caregiver hitting their shins when moving the cart.

Locking casters hold the cart firmly in place.





## SURFACE MATERIAL KYDEX 3D LAMINATE

The Guardian Delivery Cart utilizes KYDEX® 3D laminate on the top and base, which is 4x thicker than standard rigid thermofilm, giving it superior impact resistance and durability. KYDEX is certified GREENGUARD® Gold. KYDEX withstands degradation from a wide range of hospital cleaners and disinfectants.

### DURABILITY

The Guardian Delivery Cart's horizontal surfaces are thermofilmed with KYDEX® which is 4x thicker than standard thermofilm giving it superior impact strength and durability. KYDEX is certified GREENGUARD® Gold. KYDEX withstands a wide range of cleaners and disinfectants.

### INFECTION PREVENTION

The Guardian Delivery Cart has seamless horizontal and vertical panels that may aid in ease of cleaning. The large drawer has full extension drawer slides to provide easy access for cleaning.

### EASE OF USE

The pass-through drawer can be opened from front or rear of the cart. Large drawer makes delivery packs easily accessible.

Dual drop-leaf top expands caregiver workspace when needed.

Almost effortless mobility on 5" casters – 2 locking and 2 non-locking. The caster locks hold the cart firmly in place.

The curved base reduces the chances of caregivers hitting their shins while moving the cart.

### SPECIFICATIONS

Surface Width	29in (66cm)	Large Drawer Dimensions	24in W x 18in L x 13.5in D
Surface Height	36in (91cm)	Pass-through Drawer	24in W x 18in L x 5in D
Surface Depth	22in (61 cm)	Drop Leaf Extended Width	43in (1.1 m)