



HIGHLIGHTS

- **Developed for 5th wheel beds**
- **Cleanable seamless top surfaces**
- **Durable KYDEX top surface**
- **Patented height adjustment**
- **5-year gas spring warranty**
- **Manufactured in the USA**



ERGONOMIC DESIGN FEATURES

- Patent Pending base allows proper positioning with 5th wheel beds.
- Patented column facilitates smooth height adjustment.
- Activation lever is conveniently located for patient and caregiver.
- Guardian Overbed Tables are non-handed. They can be used on either side of the bed.



DESIGNED FOR 5TH WHEEL HOSPITAL BEDS

MedViron's Guardian Overbed Table is for critical care beds and stretchers with 5th wheel motorized drive systems. The Guardian Overbed Table is the only overbed table that can be properly positioned with these beds and stretchers for patient comfort and safety.

DURABILITY

MedViron's patented aluminum height-adjustment column offers superior strength, stability, and smooth height adjustment. MedViron's column has been tested for 10,000 cycles with no effect on performance.

Guardian Overbed Tables have a limited lifetime warranty. The pneumatic spring has a 5-year warranty.

INFECTION PREVENTION

A hospital's Environmental Services department is critical for infection prevention. Surfaces that are easy to clean are more likely to be cleaned. Seamless thermoformed KYDEX withstands commonly used hospital cleaners and disinfectants. The top surface is supported underneath with a steel plate which is free from obstructions and is easily wiped down.

EASE OF USE

The Guardian Overbed Table is easily positioned by a caregiver for convenient use by a patient.

The top surface can be raised/lowered with one hand. The top surface also has an eased edge for patient comfort.

The height-activation lever is conveniently located for patient and caregiver use and is used only for lowering table height.

MedViron Guardian Overbed Tables ship fully assembled in recyclable packaging.

SPECIFICATIONS

Top surface height range	33.5in (85 cm) to 45.4in (1.15 m)
Top surface dimensions	19in (48 cm) deep x 32in (81 cm) wide
Base height	.75in (1.9 cm)
Working weight capacity	25lbs (11.4 kg)
Options	Cupholders (2), Vanity mirror (1), Tablet holder (1)
Warranty	Limited Lifetime warranty, 5yr warranty on pneumatic spring