



HIGHLIGHTS

- **Cleanable seamless surfaces**
- **Durable KYDEX**
- **Spill containment top surface**
- **Robust stainless steel push handle**
- **Solid or woodgrain color options**
- **Manufactured in the USA.**



ERGONOMIC DESIGN FEATURES

- **Stainless-steel top-mounted push handle.**
- **Large open space for patient's personal belongings.**
- **Two storage drawers for caregiver supplies.**
- **Removable drawer liners to aid in cleaning.**
- **Three inch concealed casters for easy mobility.**



SURFACE MATERIAL KYDEX 3D LAMINATE

MedViron's Patient Care Cart utilizes KYDEX® 3D Laminate, which is 4x thicker than standard thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD® Gold. KYDEX withstands a wide range of cleaners and disinfectants.

DURABILITY

MedViron's Patient Care Cart is a tactical cart for caregivers. Handy and maneuverable, it's ready to accomplish tough tasks.

The robust stainless steel push handle is securely integrated with the top surface.

The Patient Care Cart is thermoformed with KYDEX, giving it superior impact strength and durability. KYDEX also withstands degradation from a wide range of cleaners and disinfectants.

INFECTION PREVENTION

A hospital's Environmental Services department is critical for infection prevention. Surfaces that are easy to clean are more likely to be cleaned. The MedViron Patient Care Cart has seamless horizontal and vertical panels which can be easily wiped down.

The drawers have 6" deep removable plastic drawer liners that may facilitate cleaning procedures.

The stainless-steel bow drawer pulls and the stainless-steel push handle can be quickly wiped down and disinfected.

EASE OF USE

Patient's belongings are often left behind if not immediately visible at departure. The Patient Care Cart has a large open space for a patient's personal items, leaving items visible and less likely to be forgotten.

Two large drawers provide ample storage and can be used to keep caregiver supplies close to the bedside.

The Patient Care Cart's spill containment area will retain 1 liter of liquid (34oz).

MedViron Patient Care Carts ship fully assembled in recyclable packaging.

SPECIFICATIONS

Cart dimensions	19.5in wide x 20.3in deep x 32.75in high (49.5 cm x 51.6 cm x 83.2 cm)
Cart height including handle	39.1 in (99.3 cm)
Cart open shelf	8.8in high (22.4 cm)
Drawer internal dimensions	15.5in wide x 14.75in deep x 5.625in high (39 cm x 37.5 cm x 14.3 cm)
Warranty	Limited lifetime warranty.