



HIGHLIGHTS

- Cleanable seamless surfaces
- Durable KYDEX®
- Impact resistant PETG tub
- Removable drawer liners
- Solid and woodgrain color options
- Manufactured in the USA

ERGONOMIC DESIGN FEATURES

- Base panel reduces the chances of caregivers hitting their shins.
- Full-extension drawer slides.
- Push handle height enables proper pushing posture.
- Patented tub design allows for easier access to patient.
- Infant tub can be inclined for medical procedures.





SURFACE MATERIAL: KYDEX® 3D LAMINATE

The Guardian Bassinet utilizes KYDEX thermoplastic on all horizontal surfaces, which is 4x thicker than standard thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD® Gold. KYDEX withstands degradation from a wide range of hospital cleaners and disinfectants.

DURABILITY

KYDEX thermoplastic, with its superior impact resistance and durability, is used on all horizontal surfaces, including the top, base panel and chart board.

The PETG plastic tub has 4 times greater impact resistance than acrylic. The tub's corners are reinforced with additional material to reduce the likelihood of corner breakage.

Guardian Bassinets have stainless steel push handles.

INFECTION PREVENTION

Seamless KYDEX surfaces withstand commonly used hospital cleaners and disinfectants.

Our stainless steel push handles eliminate the hard to get to making all of the surface area accessible and free from obstructions. This allows easy cleaning.

Our removable drawer liners have rounded corners that are easily wiped down and cleaned.

EASE OF USE

The patented tub design allows for easier access to the patient. The Infant tub can be set in an inclined position during medical procedures.

The base panel is designed to reduce the chances of caregivers hitting their shins.

A groove in the chart board keeps items from rolling off.

The bassinet's drawers have full-extension drawer slides.

SPECIFICATIONS

| | |
|---------------------|-------------------------|
| Surface Height | 30in (76 cm) |
| Handle Height | 36in (92 cm) |
| Cabinet Width | 34in (86 cm) |
| Cabinet Depth | 21.75in (53 cm) |
| Mattress Dimensions | 12.5in Wide x 26in long |
| Weight | 112lbs (51 kg) |