



**HIGHLIGHTS**

- Cleanable seamless surfaces
- Durable KYDEX
- Spill containment top surface
- Drawer-Door or 3-drawer configurations
- Solid or woodgrain color options
- Manufactured in USA

**ERGONOMIC DESIGN FEATURES**

- Removable drawer liners for ease of cleaning.
- Two stainless steel side-mounted push handles.
- Cable management grommet.
- Optional side computer access door.
- 5" twin wheel casters give the cart virtually effortless mobility.

**SPECIFICATIONS**

<b>Spill containment top</b>	20.5in wide x 21.75in deep (52 cm x 55.3 cm)
<b>Cart height</b>	36.75 in (93.3 cm)
<b>Cart width</b>	25.4in (64.5 cm)
<b>Drawer internal dimensions</b>	15.5in wide x 14.75in deep x 5.625in high (39 cm x 37.5 cm x 14.3 cm)
<b>Shelf dimensions</b>	18.5in wide x 17.75in deep (47 cm x 45 cm)
<b>Warranty</b>	Limited Lifetime warranty



### KYDEX® 3D LAMINATE

The EFM Cart utilizes KYDEX® 3D laminate, which is 4x thicker than standard thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD® Gold. KYDEX withstands a wide range of healthcare cleaners and disinfectants.

### DURABILITY

The EFM Cart's surfaces are thermofilmed with KYDEX®, which is 4x thicker than standard thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD® Gold. KYDEX withstands a wide range of cleaners and disinfectants.

Two robust stainless steel push handles provide a secure grip when moving the cart.

### INFECTION PREVENTION

A hospital's Environmental Services department is critical for infection prevention. The surfaces of the EFM Cart are completely accessible for environmental services staff to clean.

The EFM Cart has seamless horizontal and vertical panels which can be easily wiped down.

The drawers have 6" deep removable plastic drawer liners that may facilitate cleaning procedures.

Stainless steel bow drawer and door pulls are easily cleanable.

### EASE OF USE

The spill containment area of the top surface retains 20oz (590 mm) of liquid.

Drawer/Door or Three Drawer options can be made with or without key locks.

Drawer/Door cart comes with 2 plastic storage bins for instruments or supplies.

The optional side access door provides quick access to computers or power supply.

The EFM Cart has 5" TENTE® casters, 2 swivel and 2 locking, allowing the cart to be moved with little effort.

Ships fully assembled in 100% recyclable packaging.



MADE IN U.S.A.