



HIGHLIGHTS

- **Developed for 5th wheel beds**
- **Cleanable seamless top surfaces**
- **Durable KYDEX top surface**
- **Patented height adjustment**
- **5-year gas spring warranty**
- **Manufactured in the USA**

ERGONOMIC DESIGN FEATURES

The Sotto Overbed Table was designed for use with critical care beds and stretchers with 5th wheel motorized drive systems.

- Patented column facilitates smooth height adjustment.
- Activation lever is conveniently located for patient and caregiver.
- Sotto Overbed Tables are non-handed. They can be used on either side of the bed.

OPTIONS

- 32in or 39in wide surface
- Cupholders
- Vanity mirror
- Tablet holder
- Hanging shelf
- Dual Spill Containment Top (39" only)
- Two Shelf Hanging Ped (39" only)





DESIGNED FOR 5TH WHEEL HOSPITAL BEDS

MedViron's Sotto Overbed Table is specifically designed for critical care beds and stretchers with 5th wheel motorized drive systems. It is the only overbed table that can be correctly positioned with these beds and stretchers, ensuring optimal patient comfort and safety.

DURABILITY

MedViron's patented aluminum height-adjustment column provides superior strength, stability, and smooth height adjustment. This column has been independently tested for 10,000 cycles without any performance degradation.

Sotto Overbed Tables come with a limited lifetime warranty, and the pneumatic spring is covered by a 5-year warranty.

INFECTION PREVENTION

The Environmental Services department in hospitals plays a crucial role in infection prevention. Surfaces that are easy to clean are more likely to be properly maintained. Seamless thermoformed KYDEX is resistant to commonly used hospital cleaners and disinfectants. The top surface is supported underneath by an unobstructed steel plate, ensuring it can be easily wiped down.

EASE OF USE

The table can be easily positioned under a 5th-wheel bed by a caregiver for the patient's convenience.

The top surface can be raised or lowered with one hand and features an eased edge for patient comfort.

The height-activation lever, conveniently located for both patient and caregiver use, is utilized only for lowering the table height.

Tables ship fully assembled in recyclable packaging.

SPECIFICATIONS

Top surface height range	29.5in (75 cm) to 45.4in (1.15 m)
Top surface dimensions	19in (48 cm) deep x 32in (81 cm) wide, or 19in (48 cm) deep x 39in (99 cm) wide
Base height	.75in (1.9 cm)
Working weight capacity	25lbs (11.4 kg)
Options	Cupholders (2), Vanity mirror (1), Tablet holder (1), Hanging Shelf (1)
Warranty	Limited Lifetime warranty, 5-year warranty on pneumatic spring

