



MedViron 36" U-base Overbed Table



MedViron's 36" U-Base Overbed Table's ergonomic design prioritizes both patient comfort and caregiver convenience. The table's smooth, height-adjustment mechanism allows for effortless adjustment to the ideal height, reducing strain on patients and making it easy for caregivers to position the table without excessive bending or lifting. The U-shaped base makes this hospital overbed table compatible with hospital beds, patient room chairs and wheelchairs.

Options: Cupholders, Vanity Mirror or Tablet Holder, Undermount Shelf



MedViron 36" U-base Overbed Table



SPECIFICATIONS

Height Adjustment	27.5" (70 cm) to 43.75" (1.11 m)
Top Dimensions	36in wide x 19in deep (91.5cm x 48.25cm)
Base Dimensions	35in wide x 17.2in deep (89cm x 44cm)
Options	Cupholders, Vanity Mirror or Tablet Holder, Undermount Shelf
Standard Top Color Options	Honey Maple, Natural Maple, Summer Flame, Walnut, Eiderdown, Bright White
Column	Anodized aluminum
Base Color	Powder-coated steel, Arctic Silver (standard)
Casters	1.25in (32mm) black nylon twin wheel casters
Working Weight Capacity	25 lbs (11.4kg)
Weight	49lbs (22.2 kg), table only without packaging
Warranty	5 year gas spring warranty. Limited lifetime warranty for other components.