



## MedViron 36" U-base Overbed Table



MedViron's 36" I-Base Overbed Table's ergonomic design prioritizes both patient comfort and caregiver convenience. The table's height-adjustment mechanism allows for one-handed adjustment to the ideal height, reducing strain on patients and making it easy for caregivers to position the table without excessive bending or lifting. The I-shaped base makes this hospital overbed table compatible with both beds and patient room chairs.

**Options:** Cupholders, Vanity Mirror or Tablet Holder, Undermount Shelf, Drop-Leaf Extension



## MedViron 36" U-base Overbed Table



### SPECIFICATIONS

<b>Height Adjustment</b>	27.5" (70 cm) to 43.75" (1.11 m)
<b>Top Dimensions</b>	36in wide x 19in deep (91.5cm x 48.25cm)
<b>Base Dimensions</b>	35.5in wide x 18in deep (89cm x 46cm)
<b>Options</b>	Cupholders, Vanity Mirror or Tablet Holder, Undermount Shelf
<b>Standard Top Color Options</b>	Honey Maple, Natural Maple, Summer Flame, Walnut, Eiderdown, Bright White
<b>Column</b>	Anodized aluminum
<b>Base Color</b>	Powder-coated steel, Arctic Silver (standard)
<b>Casters</b>	1.25in (32mm) black nylon twin wheel casters
<b>Working Weight Capacity</b>	25 lbs (11.4kg)
<b>Weight</b>	49lbs (22.2 kg), table only without packaging
<b>Warranty</b>	5 year gas spring warranty. Limited lifetime warranty for other components.