


## MedViron Mobile Laptop Cart



### FEATURES

- Height-adjustable from 30in to 46in.
- Seamless KYDEX top surface
- Patented anodized aluminum height adjustment column
- Dual paddle height adjustment
- Curved base panel has ample room for your feet
- 6 standard colors available
- Locking casters

## MedViron Mobile Laptop Cart



### SPECIFICATIONS

<b>Height Adjustment</b>	30in to 46in (76 cm to 1.17 m)
<b>Surface Dimensions</b>	23in wide x 17in deep (58 cm x 43 cm)
<b>Base Dimensions</b>	24.5in wide x 17in deep (62 cm x 43 cm)
<b>Standard Top Color Options</b>	Honey Maple, Natural Maple, Summer Flame, Walnut, Eiderdown, Bright White
<b>Column</b>	Anodized aluminum
<b>Base Color</b>	Powder-coated steel, Arctic Silver (standard)
<b>Casters</b>	2in (50mm) twin wheel casters. Two locking, two non-locking.
<b>Working Weight Capacity</b>	25 lbs (11.4kg)
<b>Weight</b>	42lbs (19 kg) Cart only without packaging.
<b>Warranty</b>	Limited lifetime warranty.