



HIGHLIGHTS

- Cleanable seamless surfaces.
- Durable KYDEX® finish on horizontal surfaces & drawer fronts.
- The standard cart can accommodate fetal monitors that are 20" wide x 8" high. Talk to us about customization.
- Optional key locks or punch code locks.
- Solid colors or woodgrain options available.



FEATURES

MedViron's Fetal Monitor Cart's design facilitates bedside documentation, fetal monitoring, and the storage and organization of equipment and supplies. MedViron Fetal Monitor Carts have been designed to be durable and easy to clean.

MedViron Fetal Monitor Carts can be customized in size and color. The carts can be built to fit the exact requirements for the latest hardware available. They can also be manufactured to match MedViron bassinets, delivery case carts or bedside cabinets to create a complete maternity suite.

MedViron's Fetal Monitor Cart's feature an integrated pull-out keyboard tray. Single and dual height-adjustable monitor mounts can also be added.



Top Surface Dimensions	23.75" wide x 25.5" deep
Surface Height	44" high
Base Dimensions	27.5" wide x 25.4" deep
Drawer Internal Dimensions	19.5" wide x 18" deep x 6.75" high
Casters	5" diameter
Hardware	Stainless Steel Bow Pulls
Warranty	Limited Lifetime Warranty

DURABILITY

MedViron Fetal Monitor Carts are manufactured with KYDEX 3D laminate, which is 4x thicker than standard thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD Gold. KYDEX withstands degradation from a wide range of hospital cleaners and disinfectants.

INFECTION PREVENTION

The Fetal Monitor Cart has seamless horizontal and vertical panels that may aid in ease of cleaning.

The integrated pull-out keyboard tray is thermofilmed with KYDEX, making it durable and easily cleaned.

EASE OF USE

The standard cart has three large drawers for easy storage of equipment and supplies.

5" poly casters allow the cart to be rolled almost effortlessly. There are two locking casters in front and two non-locking in the rear.

The carts ship fully assembled in recyclable packaging.

STANDARD COLORS



Walnut



Summer Flame



Honey Maple



Natural Maple



Eiderdown



Bright White