



HIGHLIGHTS

- Cleanable seamless surfaces
- Durable KYDEX® finish
- Patented height adjustment
- Pass-through drawers and work surface
- Solid colors or woodgrain options available
- Manufactured in USA



ERGONOMIC DESIGN FEATURES

The Rise 2 Hospital Bassinet features height-adjustability to suit the exact needs of every user, including recovering C-section mothers. Our patented columns ensure effortless height adjustment. Each corner is equipped with independent patented height adjustment handles, enabling caregivers to easily raise or lower the bassinet's surface to the ideal level.

The innovative base panel design allows caregivers to stand closer to the work surface, protects their shins, and ensures better access and comfort while caring for newborns.

Rise 2 Bassinets can be manufactured with one or two drawers.



SURFACE MATERIAL: KYDEX® 3D Laminate

Rise 2 Bassinets are manufactured with KYDEX 3D laminate, which is 4x thicker than standard thermofilm, giving it superior impact strength and durability. KYDEX is certified GREENGUARD Gold. KYDEX withstands degradation from a wide range of hospital cleaners and disinfectants.

DURABILITY

MedViron's patented aluminum columns offer superior strength, stability, and smooth height adjustment.

The pneumatic spring has a 5-year warranty.

The PETG plastic tub has 4 times greater impact resistance than acrylic tubs.

The pass-thru storage drawer is also made with KYDEX® giving it superior durability.

INFECTION PREVENTION

A hospital's Environmental Services department is critical for infection prevention. MedViron's stainless steel push handle design makes the top surface of the bassinet completely accessible for staff to clean.

The bassinet's pass-thru drawer has a removable drawer tub with rounded inside corners for easy cleaning.

The bassinets have seamless horizontal and vertical panels which also aid in the ease of cleaning.

EASE OF USE

Patented height adjustment: Four handles, one in each corner, can be used independently to raise or lower the bassinet's height.

The pass-thru pull-out work surface can be accessed from either side of the bassinet and has a groove to keep items from rolling off.

5" casters: one swivel, two locking, one steering and locking.

Ships fully assembled in recyclable packaging.

SPECIFICATIONS

Overall Width	37.5" W
Top Surface Dimensions	36" W x 21.25" D
Surface Height	30.5" low to 40.25" high
Base Dimensions	35" W x 22" D
Drawer Internal Dimensions	22.5" W x 13.25" D x 6" H
Mattress Dimensions	26" L x 12.5" W
Hardware	Stainless Steel Bow Pulls
Warranty	Limited Lifetime Warranty