

Dura-Tough Hospital Bedside Cabinets

Durable, Easy-to-Clean Storage for Patient Rooms

MedViron's Dura-Tough Hospital Bedside Cabinet is designed for long-term performance in demanding healthcare environments. Built with durable melamine casework and a KYDEX top surface, drawer and door fronts, this medical-grade bedside cabinet resists impact, moisture, and harsh hospital cleaning agents.

The rugged construction minimizes damage from daily use, helping healthcare facilities avoid frequent repairs or costly replacements.



Features:

- KYDEX 3D Laminate tops and drawer fronts offer exceptional durability.
- Optional spill containment top or solid surface top.
- Easy to clean stainless steel drawer pulls.
- Removable drawer liners for fast, efficient cleaning.
- Optional slam lock available.
- Rolls effortlessly on twin-wheel casters.
- Limited lifetime warranty.
- Available in 8 standard colors., see our website.



customerservice@medviron.com

844-633-8476

medviron.com

Dura-Tough Bedside Cabinets



The non-porous, easy-to-clean surfaces support Environmental Services staff, while smooth-gliding drawers, ergonomic pulls, removable drawer tubs, and non-marking dual-wheel casters simplify cleaning, organization, and repositioning.

Available in three versatile configurations, Dura-Tough Hospital Bedside Cabinets are the best hospital furniture solution for facilities seeking durability, superior mobility, and easy maintenance in every patient care setting.



One Drawer/Open Shelf



Drawer/Door



Three Drawer

Specifications:

Height	30.5 inches (77.5cm)
Width	18.75 inches (47.6cm)
Depth	19 inches or (48.3cm)
Casters	37mm twin-wheel
Drawer Liner Internal Dimensions	15.375in W x 14.75in L x 5.5in D (39cm W x 37.5cm L x 14cm D)
Warranty	Limited Lifetime Warranty



customerservice@medviron.com
844-633-8476
medviron.com