

Height-Adjustable Medical Procedure Cart



Efficient, durable, and designed with caregivers in mind, the Height-Adjustable Medical Procedure Cart adapts to the unique needs of each healthcare environment. The cart is semi-customizable, allowing its width and depth to be manufactured to the specific dimensions a facility requires, making this cart one of the most flexible solutions available. From the adjustable-height work surface to the heavy-duty drawers and KYDEX surfaces, every feature is focused on improving workflow, ergonomics, and cleanability while standing up to the demands of daily hospital use.



customerservice@medviron.com

844-633-8476

medviron.com

Height-Adjustable Medical Procedure Cart

Features:

- Semi-customizable, width and depth can be modified to fit your facility's needs
- Dual height adjustment levers that can operate independently or together
- Heavy-duty interlocking drawer glide system supports up to 180lbs per drawer
- Integrated push/pull handles
- KYDEX 3D laminate exterior for long-term durability
- HDPE top surface and drawer boxes for ease of cleaning
- Ships fully assembled in 100% recyclable packaging.



Specifications:

Top Surface Height Adjustment	32.75in to 48.5in (83cm to 1.23m)
Top Surface Width as shown above	35.5in (90 cm)
Top Surface Depth as shown above	19.4in (49.3 cm)
Secondary Top Height	30.5 in (77.5 cm)
Base Width as shown above	30.3in (77cm)
Base Depth as shown above	19.4in (49.3cm)
Casters	5in (12.7cm)
Weight	196 pounds (89kg)
Warranty	Limited Lifetime Warranty



customerservice@medviron.com
844-633-8476
medviron.com