

Value Hospital Bedside Cabinets

The Value Hospital Bedside Cabinet Series delivers four versatile, budget-friendly patient room storage solutions designed for healthcare facilities seeking smart, cost-effective furniture. Each cabinet combines durable construction, ease of cleaning, and ergonomic design—perfect for keeping patient rooms organized.



Features:

- Melamine surfaces resist scratches, moisture, and everyday wear.
- Optional spill containment top or solid surface top.
- Easy to clean stainless steel drawer pulls
- Removable drawer liners for fast, efficient cleaning.
- Optional slam lock
- Rolls effortlessly on twin-wheel casters.
- Limited lifetime warranty
- Available in 8 standard colors.



customerservice@medviron.com
844-633-8476
medviron.com

Value Bedside Cabinets



Flat-pack Value Bedside Cabinet

All melamine, t-mold melamine top,
open case with one fixed shelf.
Ships knock-down, assembly required.

Fully-assembled Value Bedside Cabinets

Fully assembled Value configurations are manufactured
with a melamine case with plastic drawer glides as standard:

Single drawer/open fixed shelf.
One drawer and one door.
Three drawers.



Specifications:

Height	30.5 inches (77.5cm)
Width	18.75 inches (47.6cm)
Depth	19 inches or (48.3cm)
Casters	37mm twin-wheel
Drawer Liner Internal Dimensions	15.375in W x 14.75in L x 5.5in D (39cm W x 37.5cm L x 14cm D)
Warranty	Limited Lifetime Warranty



customerservice@medviron.com
844-633-8476
medviron.com