

Martello Two-Drawer Hospital Bedside Cabinet



Features

The Martello Two-Drawer Hospital Bedside Cabinet was designed for both patients and staff, featuring a large lower drawer for backpacks, purses and personal items, and a front-mounted USB charging port. Features a solid surface top for durability and everyday performance, and a removable top drawer tub for easy cleaning. Rolls effortlessly on hidden twin-wheel casters.

- Solid surface top
- USB-A & USB-C charging ports
- Concealed twin-wheel casters
- Removable top drawer plastic liner for easy cleaning
- Large bottom drawer for backpack and purse storage
- KYDEX thermoplastic exterior finish
- Limited lifetime warranty



customerservice@medviron.com
844-633-8476
medviron.com

Martello Two-Drawer Hospital Bedside Cabinet



Specifications

Height: 30 in (nominal)

Width: 18.75 in

Depth: 19 in

Weight: 67 lb

Casters: 37mm Twin-wheel

Bottom Drawer ID: 16in W x 15.3in D x 9.4in H

USB Charging: Front-mounted USB-A and USB-C charging ports with 20W charging output and a standard 6 ft, 120V power cord.

Large bottom drawer for backpacks or purse storage



Integrated USB charging ports



Removable drawer tub for ease of cleaning



customerservice@medviron.com

844-633-8476

medviron.com